


Pack/unpack sample

This feature is present in ADCI v1.5 or newer

In the event a [sample](#) must be moved from one computer to another, multiple components of the sample must be transferred to the destination computer. Sample packing/unpacking provides functionality to simplify this process.

Pack (export) sample



Unpack (import) sample

Sample Packing/Unpacking

Sample Information

Sample ID:

Description:

This is an example sample used to demonstrate ADCIDE functionality. It should not be used for any real tests.

It includes only 8 images allowing it to be processed in a short period of time.

Our sincere gratitude to Dr. Ruth Wilkins, Radiobiology and Protection Division, Health Canada, for providing these images.

Contour files: found. Metaphase images: found.

Unpacking

Unpack Sample and Save As File:

Specify File

☒ Unpack Contour Files and Copy to ADCI Data Folder

Unpack metaphase images and copy to specified folder:

Specify Folder

Progress

Ready to unpack sample

0%

Start Unpacking